

COMPONENT SHEET

K-1615 COMBINATION LAUNDRY WASHROOM CONTROL

COMPONENTS (not Bill of Materials):

P/T Alkalinity

1 x 5013	Instruction
1 x 9198G	Sample Tube, Graduated (25 mL) w/ cap & green dot, plastic
1 x R-0638G-C	Phenolphthalein Indicator, 2 oz w/ green cap, DB
1 x R-0645-C	Total Alkalinity Indicator, 2 oz, DB
1 x R-0687G-C	Sulfuric Acid .12N, 2 oz w/ green cap, DB
1 x R-0736G-C	Sulfuric Acid .6N, 2 oz w/ green cap, DB

Chlorine (Bleach)

1 x 4026	Dipper Spoon, 2 g, plastic, white
1 x 4078	Pipet (eye dropper), Graduated (3 mL w/ 0.5 mL div), plastic
1 x 5004	Instruction
1 x 9009	Pipet (eye dropper), Calibrated (0.5 & 1.0 mL) w/ yellow cap, plastic
1 x 9198Y	Sample Tube, Graduated (25 mL) w/ cap & yellow dot, plastic
1 x R-0664-C	Bleach Reagent #1, 2 oz
1 x R-0665S-II	Bleach Reagent #2 (crystals), 50 g
1 x R-0666-C	Bleach Reagent #3, 2 oz, DB

Total Hardness

1 x 5202	Instruction
1 x 9198B	Sample Tube, Graduated (25 mL) w/ cap & blue dot, plastic
1 x R-0619B-C	Hardness Buffer, 2 oz w/ blue cap, DB
1 x R-0620B-I	Hardness Indicator Powder, 10 g w/ blue dot
1 x R-0683-C	Hardness Reagent (1 drop = 10 ppm hardness as CaCO ₃), 2 oz, DB

pH

1 x 5425	Instruction
1 x 9017	Test Cell, Calibrated (5 mL) w/ cap, plastic
1 x R-1003U-C-DB	Long Range Indicator, 2 oz, DB

Misc.

1 x 5063	Instruction (Component Sheet)
----------	-------------------------------

Instr. #5063

TO ORDER REPLACEMENT PARTS AND REAGENTS CALL TOLL-FREE
800-TEST KIT (800-837-8548).

**CAREFULLY READ AND FOLLOW PRECAUTIONS ON REAGENT LABELS.
KEEP REAGENTS AWAY FROM CHILDREN.**

NOTE: When dispensing reagents from dropper bottles, **always** hold bottle in a vertical position.



31 Loveton Circle, Sparks, MD 21152 USA
800-TEST KIT (837-8548) • 410-472-4340

061918